

# Biomass Compositional Analysis

## New Tools and Methods Supporting Biomass Utilization

May 1-2, 2003

Bonnie R. Hames  
National Renewable Energy Laboratory

# Biomass Analysis at NREL



Biomass Analysis Experts: more than 20 years of experience  
Experienced in the analysis of a wide variety of biomass types  
Standard Methods published through ASTM E48 Biotechnology  
Standard Reference Materials available through NIST.

# Importance of Biomass Analysis in Biomass Utilization

Natural heterogeneity of biomass

- Complex polymer matrix
- Inherent property of biomass
- Difficult to control in feedstock
  - Function of many variables
  - Feedstocks often residues of other process


Significant challenge for emerging industries

- Complexity of biomass underestimated
- Traditional methods are expensive and slow
- New tools and methods needed

# Advantages of Rapid Biomass Analysis

- Useful in industrial applications
  - Demonstrated technology
  - High impact / low risk
- Faster
  - Minutes instead of days
  - Minimal sample preparation
- Cheaper
  - \$10 - \$20 per sample
  - Compared to \$800-\$2,000 for wet analysis
- Better
  - Calibrated using best methods
  - Less operator dependent

# Applications of Rapid Biomass Analysis




Real-time  
monitoring

Carbohydrates  
Lignin  
Protein  
Moisture  
Cell Mass

Ethanol  
Other Products  
Heating Value  
Ethanol Production Capacity

# Integrated with Other Projects



# Rapid Biomass Analysis Method Development Essentials

## 1. Calibration Samples

- Determines the range of new method
- Robust model requires approximately 100 samples
- Should resemble samples to be analyzed
- Account for variability (region, season, process, severity etc.)
- Easier to obtain in research environments
- 

## 2. Chemical Characterization

- Determines precision and accuracy of new method
- Requires appropriate analytical methods
- Most expensive part of method development

# Rapid Biomass Analysis

## Method Development Essentials

### 3. Rapid Technique

- Determines speed and cost of new method
- Should be robust and reproducible in collection environment
- Must be sensitive to compositional changes

### 4. Multivariate Analysis

- Translates spectroscopic data into compositional information


### 5. QA/QC

- Guides use of new method
- Calibration checks (include use of well characterized standard reference materials)
- Sample screening (use of outlier flags, e.g.  $GH=3.0-2.5$  std dev 99.5% population)



# Improving Calibration Information

## Corn Stover Composition



# Wet Chemical Methods for Corn Stover Analysis


## Methods in Use

- Carbohydrates
  - Glucan
  - Xylan
  - Arabinan
  - Mannan
  - Galactan
- Lignin
  - Acid Insoluble
  - Acid Soluble
- Protein (N%)
- Soil
- Structural Inorganics
- HMF / Furfural
- Acetyl groups


## Methods being Developed or Improved

- Sucrose
- Uronic Acids / Pectins
- Amino Acids / Protein
- Nitrates / Nitrites
- Chlorophyll
- Protein Degradation Products
- Carbohydrate Degradation Products
- Physical properties


# Corn Stover Feedstock NIR/PLS Method




# Pretreated Corn Stover : Dry solids




# Pretreated Corn Stover: Slurry Method



# Predicted vs. Measured Pretreatment: liquors1.eqa



- Sugars
  - Monomers and Oligomers
- Glycerol
- pH
- Acetic Acid
- HMF and Furfural




- Transmission spectroscopy
- Helma, flow-through quartz cuvette cell (pathlength of 0.5mm)
- Water used as reference
- Temperature Controlled at 60°C

# Investigating Different Spectroscopic Techniques

Dry solid method

- Reflectance spectroscopy
- Slurry method
- Transflectance spectroscopy
- Liquid method
- Transmission spectroscopy



# Improving Research Providing New Levels of Information

- Rapid Analysis Methods allow experiments that require analysis of hundreds of samples
  - Not possible using traditional methods
    - Too expensive by wet chemical methods
      - \$2,000/ sample x 1,000 samples = \$2,000,000
    - NIR/PLS
      - \$20/sample x 1,000 samples = \$20,000
  - Too slow by wet chemical methods
    - Years to process 1,000 samples
  - NIR/PLS
    - 5-7 days to analyze 1,000 samples



# Plans for Improving Wet Chemical Methods

- Expand slate of analyses
  - New constituents
  - New feedstocks
  - New processes
  - Physical properties measurements
- Improve current methods
  - Improve precision and accuracy
  - Improve total mass closures
  - Expand QA/QC methods
- Develop and validate feedstock-specific methods
  - Publish Validated methods through ASTM
  - Post on NREL website

# Plans for Future Rapid Analysis Methods Development

- Demonstration in NREL Pilot Plant
  - On-line and at-line methods
  - Calibration transfer protocols
- Improve Field methods
  - Feedstock assessment
  - Genetic screening
- Expand slate of Rapid Analysis Methods
  - Different pretreatment chemistries
  - Corn Fiber, Corn Bran, DDG
  - SSF residues and fermentation samples
  - Field methods
  - Elemental Analysis and Higher Heating Value of process residues

# Critical Issues

- High Cost
  - Cost of RA method development around \$300,000
  - Methods are feedstock and process specific
  - Many methods needed to monitor entire biomass conversion process
- Resource issues
  - Current pace and scope of these projects cannot provide rapid analysis methods for industrial scale on-line applications in the near future

# Biomass Rapid Analysis Network BRAN

- Consortium for Rapid Biomass Analysis Method Development
  - Fast-track rapid analysis method development
  - Expand calibrations for general applications
  - Share costs
    - Sample collection
    - Wet chemical analysis
    - Spectroscopic analysis
  - Consensus method validation
  - Creating and training future workforce
  - Communication network

# Biomass Analysis

## Acknowledgements

- Wet Chemical Methods
  - David Johnson
  - Ray Ruiz
  - Chris Scarlata
  - Amie Sluiter
  - Justin Sluiter
  - David Templeton
- Rapid Analysis Methods
  - Ray Ruiz
  - Amie Sluiter
  - David Templeton
- Industrial Partnerships / BRAN
  - John Ashworth